







1	2	3	4
REV.	DR:	DATE	CRAFTWORK



细节 A
比例 2.000




Notes

1. 标⊕号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和锐角等
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 表格尺寸参照3D图

STATE: Release

Select			V	UNITS	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA		
	± A	± B	± C	MM				
				MAT'L			PART NO. (INTENDED USE)	TITLE:
	0-6	0.02	0.05	0.10				PCR20-C-96-HS-CA
	>6-30	0.05	0.10	0.15			FINISH	DWG NO.:
	>30-60	0.10	0.15	0.20				
>60-120	0.15	0.20	0.25	Q'TY	CHKD:			
>120-180	0.20	0.25	0.30		DR:			
>180-300	0.25	0.30	0.35					



SCALE	SHEET	REV.
1:1	1/1	A

A4