




REV.	DR:	DATE	NOTE	CRAFTWORK
------	-----	------	------	-----------




截面 A-A
比例 1.500



Notes

1. 标*号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 未标注尺寸参照3D档

STATE: Release

Select	UNITS			NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA						
	MM										
	± A	± B	± C					PART NO. (INTENDED USE)	TITLE: RH50-H-E/C		
	0-6	0.02	0.05 0.10								
	>6-30	0.05	0.10 0.15					APPD:	DWG NO.:		
	>30-60	0.10	0.15 0.20								
>60-120	0.15	0.20 0.25	CHKD:		SCALE	SHEET	REV.				
>120-180	0.20	0.25 0.30			DR:	1:1	1/1	A			
>180-300	0.25	0.30 0.35	Q'TY								