

A

A

B

B

C


C

D



STATE: Release

E

A4



截面 A-A

**Notes**

1. 标⌀号为重点管径尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 未标注尺寸参照3D图

Select				UNITS MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA								
	± A	± B	± C	MAT'L PP	PART NO. (INTENDED USE)	TITLE: GT125ul								
0-6	0.02	0.05	0.10	FINISH	APPD:	DWG NO.:								
>6-30	0.05	0.10	0.15	Q'TY	CHKD:	<table border="1" style="display: inline-table; border-collapse: collapse;"> <tr> <td style="text-align:center;">SCALE</td> <td style="text-align:center;">SHEET</td> <td style="text-align:center;">REV.</td> </tr> <tr> <td style="text-align:center;">1:1</td> <td style="text-align:center;">1/1</td> <td style="text-align:center;">A</td> </tr> </table>			SCALE	SHEET	REV.	1:1	1/1	A
SCALE	SHEET	REV.												
1:1	1/1	A												
>30-60	0.10	0.15	0.20		DR:									
>60-120	0.15	0.20	0.25											
>120-180	0.20	0.25	0.30											
>180-300	0.25	0.30	0.35											