



Notes

1. 标Φ为直径管径尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 未标注尺寸参照3D图

STATE: Release

Select	UNITS				NAME (INTENDED USE)	浙江友莱生物科技有限公司				
		v				MM	YOU LAI BIOTECHNOLOGY CO.,LTD.			
		± A	± B			± C	PART NO. (INTENDED USE)	ZHE JIANG, CHINA		
	0-6	0.02	0.05			0.10		TITLE:	T1000吸头	
	>6-30	0.05	0.10			0.15	APPD:		DWG NO.:	
	>30-60	0.10	0.15			0.20	FINISH	CHKD:		
>60-120	0.15	0.20	0.25	Q'TY	DR:	SCALE	SHEET	REV.		
>120-180	0.20	0.25	0.30	I		1:1	1/1	A		
>180-300	0.25	0.30	0.35							