


- Notes**
1. 标*号为重要管控尺寸
 2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕等
 3. 产品要求洁净无油 生产场所 10万级净化车间
 4. 未标注尺寸参照3D档

STATE: Release

Select	v			UNITS	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA			
	± A	± B	± C	MM					
	0-6	0.02	0.05	0.10			MAT'L		
	>6-30	0.05	0.10	0.15			PP		
	>30-60	0.10	0.15	0.20			FINISH		
	>60-120	0.15	0.20	0.25			APPD:		
>120-180	0.20	0.25	0.30	O'TY	CHKD:	PART NO. (INTENDED USE)	TITLE:		
>180-300	0.25	0.30	0.35	DR:	DWG NO.:	RR200-B-E	SCALE	SHEET	REV.
							1:1	1/1	A



A4