


截面 A-A  
比例 1.500




比例 1.500

**Notes**

1. 标\*号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和磕伤
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 未标注尺寸参照3D图

STATE: Release

Select	UNITS			MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA								
	± A	± B	± C			PART NO. (INTENDED USE)	TITLE: RT50-B-E/C							
	0-6	0.02	0.05							0.10				
	>6-30	0.05	0.10							0.15	APPD:	DWG NO.:		
	>30-60	0.10	0.15							0.20				
	>60-120	0.15	0.20							0.25				
>120-180	0.20	0.25	0.30	CHKD:										
>180-300	0.25	0.30	0.35		DR:									
				O'TY		I								
				UNITS	SCALE SHEET REV.									
				MM	1:1	1/1	A							



A4