

1

2

3

4


REV.

DR:

DATE

NOTE

CRAFTWORK



$\varnothing 18.20$   
 $\varnothing 15.30$

15.00

10ml

9

8

7

6

5

4

3

2

1

$\varnothing 11.21$

165.00

118.00

$\varnothing 2.00$

$\varnothing 2.90$

1. 未注尺寸公差: X $\pm$ 0.3; X.0 $\pm$ 0.1; X.00分三种: (0<X $\leq$ 20,  $\pm$ 0.005; 20<X $\leq$ 150,  $\pm$ 0.01; X>150,  $\pm$ 0.015); X.000 $\pm$ 0.005, 轴类用大尺寸, 孔类用小尺寸。  
2. 相同形状标识: A.B.C.D.....  
3. 未标注的螺钉深度等于螺丝直径的两倍。  
4. 未标注螺丝深度等于螺丝直径两倍。  
5. 料直径垂直度: 0<长, 宽 $\leq$ 100 0.01;  
6. 外形尺寸公差: 0.00 $^-$ 0.01; 孔尺寸公差: +0.01 $^-$ 0.00; 拖位尺寸公差: 托槽+0.01 $^-$ +0.03; 托0.00 $^-$ 0.02

UNITS

MM

NAME (INTENDED USE)

MAT'L

PP

PART NO. (INTENDED USE)

TITLE:

10ml小口吸头

APPD:


DWG NO.:

10ML-X10-2

CHKD:

DR:

08-Aug-22



SCALE

1.000

SHEET

1/1

REV.

A

1

2

3

4

A

A

B

B

C

C

D

D

E

E

A4