

A

A



B

B

C

C



D



◎	0.8	B
---	-----	---

∅	0.50	+0.03	-0.03
---	------	-------	-------

截面 A-A  
比例 1.500



**Notes**

1. 标\*号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕等
3. 产品要求洁净免洗 生产场所 10万级净化车间
4. 未标注尺寸参照3D图

STATE: Release

		V		UNITS MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA		
Select	± A	± B	± C	MAT'L	PART NO. (INTENDED USE)	TITLE:		
0-6	0.02	0.05	0.10	PP	APPD:	T200吸头		
>6-30	0.05	0.10	0.15	FINISH	CHKD:	DWG NO.:		
>30-60	0.10	0.15	0.20	Q'TY	DR:	SCALE	SHEET	REV.
>60-120	0.15	0.20	0.25	1	1:1	1/1	A	A4
>120-180	0.20	0.25	0.30					
>180-300	0.25	0.30	0.35					